# Wirobond® 280

# Wirobond® C

# Wirobond® easy

# Wirobond® LFC

# Wirobond® SG

	Wirobond® 280 Seite/page	Wirobond® C Seite/page	Wirobond® easy Seite/page	Wirobond® LFC Seite/page	Wirobond® SG Seite/page
<b>de</b> Gebrauchsanweisung	2	30	58	86	114
<b>en</b> Instruction for use	3	31	59	87	115
<b>fr</b> Notice d'utilisation	4	32	60	88	116
<b>es</b> Instrucciones de uso	5	33	61	89	117
<b>it</b> Istruzioni per l'uso	6	34	62	90	118
<b>nl</b> Gebruiksaanwijzing	7	35	63	91	119
<b>pt</b> Instruções de utilização	8	36	64	92	120
<b>ru</b> Инструкция по применению	9	37	65	93	121
<b>ar</b> دليل الاستخدام	10	38	66	94	122
<b>bg</b> Инструкция за употреба	11	39	67	95	123
<b>cs</b> Návod k použití	12	40	68	96	124
<b>da</b> Brugsanvisning	13	41	69	97	125
<b>el</b> Οδηγίες χρήσης	14	42	70	98	126
<b>et</b> Kasutusjuhend	15	43	71	99	127
<b>hr</b> Upute za uporabu	16	44	72	100	128
<b>hu</b> Használati utasítás	17	45	73	101	129
<b>ja</b> 使用説明書	18	46	74	102	130
<b>ko</b> 사용 설명서	19	47	75	103	131
<b>lt</b> Naudojimo instrukcija	20	48	76	104	132
<b>lv</b> Lietošanas pamācība	21	49	77	105	133
<b>pl</b> Instrukcja używania	22	50	78	106	134
<b>ro</b> Instrucțiuni de utilizare	23	51	79	107	135
<b>sk</b> Návod na použitie	24	52	80	108	136
<b>sl</b> Navodilo za uporabo	25	53	81	109	137
<b>sv</b> Bruksanvisning	26	54	82	110	138
<b>th</b> คำแนะนำการใช้งาน	27	55	83	111	139
<b>tr</b> Kullanma talimatı	28	56	84	112	140
<b>zh</b> 使用说明	29	57	85	113	141

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
## Instruction for use

Dental Co-based metal-ceramic alloy, Type 4. Wirobond® 280 is available as cylinders. Wirobond® 280 complies with ISO 22674 and ISO 9693. REF 50134 – 1000 g; REF 50135 – 250 g; REF 50159 – 1200 g; REF 50136 – 24 g sample	
<b>Alloy characteristics</b>	
According to ISO 22674 free of nickel, cadmium, beryllium and lead	
Type (accord. to ISO 22674)	4
Preheating temperature	°C 900–1000
Solidus, liquidus temperature	°C 1355, 1430
Density	g/cm <sup>3</sup> 8.6
Casting temperature	°C 1500
Young's modulus	GPa 221 / *214
Proof strength (R <sub>p0.2</sub> )	MPa 475 / *500
Elongation after fracture (A <sub>5</sub> )	% 9 / *5
Vickers hardness	(HV10) 280 / *390
Coefficient of thermal expansion (CTE) 25 – 500 °C, 10 <sup>-6</sup> K <sup>-1</sup> (cast / * after ceramic firings)	14.3 / *14.0
BEGO color code	8
Investment material: phosphat bonded, e. g. Bellavest SH (REF 54252)	
Crucible material: ceramic	
Melting powder: Wiromelt (REF 52526)	
Veneering ceramic Ceramic with suitable CTE, e. g.: VITA VMK Master	
Oxidation firing is not recommended, but these parameters can be used: 900 °C/5 min/vac	
Highest recommended firing temperature: 980 °C	
Heating rate recommended max. 55 °C/min	
Flux: e. g. Minoxid (REF 52530)	
Brazing material before firing: Wirobond®-Lot (REF 52622)	
Laser wire: Wiroweld (REF 50003; 50005)	

**Intended Use:** Cobalt-based metal-ceramic alloys are intended for casting of dental restorations such as crowns, bridges as well as metal-ceramic restorations.

**Indication:** Cobalt-based metal-ceramic alloys are intended to treat the condition of missing hard tissue (teeth).

**Contraindications:** No contraindications are known. However, unwanted biological reactions such as allergies to contents of the alloy or electrochemically based reactions may very rarely occur. In case of known incompatibilities and allergies to contents of the metallic material it should not be used.

**Clinical benefit:** Artificial replacement of hard tissue (teeth), to restore masticatory functionality (aesthetic and function).

**Warnings:** Metal dust is harmful to your health. Grinding and blasting should be done under a suitable vacuum. Respiratory protection of type FFP3-EN149 is recommended!

**Precautions:** Electrochemically induced sensitivity may very rarely occur when the restoration makes occlusal or proximal contact with a different alloy. This product has not been evaluated for safety and compatibility in the MR environment. It has not been tested for heating, migration, or image artifact in the MR environment. Safety in the MR environment is unknown. Scanning a patient with this device may result in patient injury.

**Patient group:** Objects made of the alloy can be used regardless of patient age. The alloy should not be used in cases of known incompatibilities or known allergies to alloy components.

**Adverse reactions:** No adverse reactions are known. Nevertheless, rare cases of individual reactions to single components can never be excluded completely. In this case, the application should not be continued.

**Wax up:** Minimum wall thickness after grinding: 0.3 mm. Avoid sharp edges and corners. Frameworks for veneering should be anatomically reduced. Connectors should be modeled as wide and high as possible (height: min. 3.5 mm, width: min. 2.5 mm). In case of bruxism stronger modellation is required. Use wax or plastic hollow sticks. Do not taper the spruing.

**Investing:** Use only phosphate bonded investment material for crowns and bridges.

**Melting / casting:** Do not overheat alloy. Use only clean ceramic crucibles and one crucible per alloy. To enable an exact batch traceability, cast new metal only. If applicable use melting powder. Follow the casting device manufacturer's instructions for parameters and casting procedures. After casting, cool the mould slowly.

**Grinding:** Use tungsten carbide burs.

**Polishing:** Polishing may be eased by blasting the surface with lead free soda glass (e.g. Perlablast® micro). Afterwards polish with suitable polishing paste using rubber polishers and brushes. Clean surface thoroughly by steam cleaning or boiling in distilled water.

**Ceramic veneering:** Use veneering ceramics with suitable CTE (ISO 9693). Follow the ceramic manufacturer's instructions for use. Remove the oxide by blasting (250 µm / 3-4 bar; e. g. with Korox 250). Clean surface thoroughly by steam cleaning or boiling in distilled water. Do not touch surfaces afterwards with hands. Use artery clamps or similar devices. Support the frameworks adequately during firing cycles.

**Soldering / brazing:** Fix the parts with soldering investment material (e. g. Bellatherm®). The prepared gap shall not exceed 0.2 mm with parallel walls. Use a suitable BEGO flux material. The flux residues and oxides must etched off. Clean surface thoroughly by steam cleaning or boiling in distilled water.

**Laser welding:** If possible use X-seam and filler material. Follow manufacturer's instructions for use and hazard notes of the laser welding devices! Clean surface thoroughly by steam cleaning or boiling in distilled water.

**Limit of Liability:** Except where prohibited by law, BEGO Bremer Goldschlägerei Wilh. Herbst GmbH & Co. KG will not be liable for any loss or damage arising from this product, whether direct, indirect, special, incidental or consequential, regardless of the theory asserted, including warranty, contract, negligence or strict liability.


**Storage conditions:** No special storage conditions are known.

**Warranty:** Whether given verbally, in writing or by practical instructions, our recommendations for use are based upon our own experience and trials and can be considered as standard values. Our products are subject to a constant further development; therefore, alterations in construction and composition are reserved. Any serious incident that has occurred in relation to the products should be reported to BEGO Bremer Goldschlägerei Wilh. Herbst GmbH & Co. KG and the competent authority.

**Instructions for disposal:** Disposal procedures as follows:

**Device:** The assignment of a waste key identification number as per the European Waste Catalogue Ordinance (AVV) must be carried out in consultation with the regional waste disposal contractor. Do not dispose of with household waste.

**Packaging:** Packaging must be fully emptied and properly disposed of in compliance with statutory regulations. Packaging that is not fully emptied must be disposed in coordination with the regional waste disposal contractor.

Consult instructions for use

Caution

Use-by-date

Batch code

Contains hazardous substance

Non-sterile

Rx only  
Only for technical  
personnel!

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