

MT3

Wet trimmer

The MT3 wet trimmer is impressive due to its exceptionally powerful motor. With the MT3 Renfert provides high-performance trimming at a low price.

Advantages

- Quick change of the trimmer table angle using the practical tipping mechanism (90° / 98°).
- Safety: when the front door is opened a safety switch stops the motor and water.
- Uniform coverage by the optimally positioned spray tube prevents clogging of the trimmer disc.

Details

- Optimum view of the model due to a 10° inclination of the unit.
- Easy guidance of the plaster model with the aid of slotted angle lines on the trimmer table.



tool in-
cluded



easy use



Uncompromisingly powerful – let it convince you:
www.renfert.com/mt3

Technical data

Permissible mains voltage	220–240 V 100–120 V
Permissible mains frequency	50/60 Hz
Power consumption	1300 W (230 V) 1325 W (120 V) // 1.74 hp (230 V) 1.78 hp (120 V)
Motor performance (P1)	1300 W // 1.74 hp 1.78 hp
Motor performance (P2)	900 W // 1.21 hp 1.34 hp
Rotational speed	2880–3400 rpm
Water consumption max.	7 l/min // 0.25 cfm
Water pressure min. / max.	1–5 bar // 14.5–72.5 psi
Ø (Connection piece for drain hose)	36 mm // 1.42 "
Ø (Trimmer disc)	234 mm // 9.2 "
Dimensions (W x H x D)	305 x 330 x 410 mm // 12.0 x 13.0 x 16.2"
Weight (without trimmer disc)	~13,2 kg // ~29 lbs

Ordering information

MT3 incl. Klettfix trimmer disc, 220–240 V	No. 18080000
MT3 incl. Marathon trimmer disc, 220–240 V	No. 18080500
MT3 incl. Klettfix trimmer disc, 100–120 V	No. 18081000
MT3 incl. Marathon trimmer disc, 100–120 V	No. 18081500
MT3 100–120 V 50/60Hz JP + Klettfix trimmer disc	No. 18082000
ORTHO guide	Page 62
Trimmer discs	Page 63