



KENDA ALL CERAMIC

2-STEP DIAMOND POLISHING SYSTEM

- Flexible polisher for all ceramic restorations
- Adapts to tooth shape
- Prevents abrasion of natural opposing teeth
- Reduces the risk of fracture and chipping



SET REF 0200.006A



19

REF	reorder	0215	0216	0218	0219	0225	0226	0228	0229
ISO	1 : 868 204 ... 2 : 838 204 ...	030 523 060	243 523 030	243 523 050	303 523 100	030 503 060	243 503 030	243 503 050	303 503 100
	Ø x L mm	6,5 x 10	3 x 7	5 x 10	10 x 1	6,5 x 10	3 x 7	5 x 10	10 x 1
μ	grit	color	rec.	max.	μ	grit	color	rec.	max.
	medium	medium bordeaux	7500 - 10 000 RPM	10 000 RPM		extra-fine	light bordeaux	5000 - 7 500 RPM	10 000 RPM
Refill Box of:	N Pressure:	Autoclavable:	Shank:	Medical device:					
6, 12, 25, 100	0,3 - 2N	135°C, 275°F	RA, stainless steel	Class IIa					

