

Cement selection made easy.

	RelyX™ Ultimate Adhesive Resin Cement	RelyX™ Unicem RelyX™ Unicem 2 Self-Adhesive Resin Cement	RelyX™ Ketac Cem Plus in Clicker	RelyX™ Veneer Cement
Metal/Metal based				
Inlays/Onlays	+	++	+	-
Crowns/Bridges	+	++	++	-
Endodontic posts	+	++	++	-
Maryland bridges	++	+	-	-
On implant abutments	+	++	++	-
Glass ceramics (incl. Li. Disilicate) (e.g. e.max®, VITA Mark II, IPS Empress 2)				
Inlays/Onlays/Table tops	++	+	-	-
Crowns/Bridges	++	++	-	-
Veneers	+	-	-	++
On implant abutments	+	+	-	-
Oxide ceramics (e.g. Lava™ Plus, Brux Zir®, Procera®)				
Inlays/Onlays	++	++	+	-
Crowns/Bridges	+	++	+	-
Endodontic posts	+	++	+	-
Maryland Bridges	++	+	-	-
On implant abutments	+	++	++	-
Resin nano ceramics (e.g. Lava™ Ultimate CAD/CAM Restorative)				
Inlays/Onlays	++	-	-	-
Veneers	+	-	-	++
Resin composites				
Inlays/Onlays	++	+	-	-
Crowns	++	++	-	-
Endodontic posts	+	++	-	-
Veneers	+	-	-	++
On implant abutments	++	++	-	-

++ Highly recommended* + Recommended - Not indicated

* Either better performance for this indication or easier handling with equal performance.

Note: The cement selection recommendations in this guide list the key indications for a particular cement, and are intended as a general guide only. 3M RelyX cements are highly versatile products that can be used for multiple indications. Your selection is dependent on individual patient circumstances and your preferred procedural techniques.

3M Oral Care

3M United Kingdom PLC
Charnwood Campus
10 Bakewell Road
Loughborough
LE11 5RB

0800 626 578

3M Ireland Ltd

The Iveagh Building
The Park, Carrickmines
Dublin D18
Ireland

+353 (0)1280 3555

Learn more at
3M.co.uk/Dental

3M, Lava, RelyX and Scotchbond are trademarks of the 3M company.
Please recycle. © 3M 2018. All rights reserved.J348521.