

# Get a great start to intraoral scanning with 3Shape TRIOS

## Technical Specifications



### TRIOS 6\*

*Revolutionary.  
TRIOS Redefined.*



### TRIOS 5\*

*Intelligent.  
Acclaimed technology.*



### TRIOS 3\*

*Proven.  
The trusted choice.*



### TRIOS Core\*

*Core.  
Essential functionality.*

Dimensions, HWL (cm)	3.8 x 3.9 x 26.6	3.8 x 3.7 x 26.6	4.8 x 4.0 x 27.3	4.8 x 4.0 x 27.3
Total weight, gram (incl. battery and tip)	308	299	340	340
Images captured by the scanner per second	2400	2400	1875	1400 (at least)
Proven clinical accuracy	✓	✓	✓	✓
Scanner and closed tip cleared for use in USA (FDA 510(K))	✓	✓	X	X
Calibration free	✓	✓	X	X
Wireless scanning	✓	✓	X	X
TRIOS Share compatible	✓	✓	X	X
ScanAssist to ensure continuous alignment of scan data	✓	✓	X	X
AI Scan 2.0 to easily scan challenging parts of the jaw	✓	✓	✓	✓
High-definition scanning	✓	X	X	X
Realistic colours and shade measurement	✓	✓	✓	✓
TRIOS Patient Specific Motion to ensure fit of restoration design	✓	✓	✓	✓
Remote Control built into the scanner	Advanced	Advanced	Simple	Simple
Compatible with MOVE Pro	✓	✓	X	X
Compatible with MOVE+	X	✓	✓	X

\* The availability of the various TRIOS models depends on market registrations. Please contact your local 3Shape sales representative regarding availability.

# Discover the 3Shape TRIOS Software and Apps

## Software and Apps



### TRIOS 6\*

*Revolutionary.*  
Acclaimed technology.

### TRIOS 5\*

*Intelligent.*  
Acclaimed technology.

### TRIOS 3\*

*Proven.*  
The trusted choice.

### TRIOS Core\*

*Core.*  
Essential functionality.

Seamless workflows with 25.000+ labs and 100+ partners via 3Shape Unite



PMS integration via 3Shape Unite



Patient Engagement apps



TRIOS Dx Plus for AI assistive detection of most common dental pathologies



TRIOS Dx Standard to measure, compare, inspect over time



Model Builder Express for print-ready model design



Invisalign connection



Implant Planner, Ortho Planner

Add-on Option

Add-on Option

Add-on Option

Add-on Option

Studio Apps for in-house design, planning and production

Add-on Option

Add-on Option

Add-on Option

Add-on Option

\* The availability of the various TRIOS models depends on market registrations. Please contact your local 3Shape sales representative regarding availability.

From TRIOS to patient outcome. Unite everything in between.

## A great scan is just the beginning

3Shape TRIOS gives you accurate scans and an unparalleled scanning experience. 3Shape Unite is the workflow engine that unlocks a world of possibilities with your scan. From digital impression to treatment conclusion, it unites everything in between.



Treatments  
made simple



Dental apps  
made accessible



Patient cases  
made available

